

## Chromium - Hexavalent and total Chromium

**Method: DTPA in acetate buffer, pH 6.2**

**Function: Differential Pulse Adsorptive Stripping Voltammetry (DPS/a)**

|                         |               |
|-------------------------|---------------|
| Start Potential (mV)    | -1000         |
| End Potential (mV)      | -1400         |
| Current range           | 20,48 $\mu$ A |
| Scan Speed (mV/s)       | 10            |
| Deposition time (s)     | 30            |
| Deposition Pot. (mV)    | -900          |
| Number of cycles        | 3             |
| Delay before sweep (s)  | 5             |
| Purge and stir time (s) | 300           |
| Stirring speed (rpm)    | 300           |
| Drop Size (a.u.)        | 60            |

### Chromium (VI) concentrate standard solution (1 g/l)

Dissolve 2.828g of pure and dried  $K_2Cr_2O_7$  in 1 l of 1%  $HNO_3$ , in a volumetric flask. ( $MM_{K_2Cr_2O_7} = 294.19$ ;  $MM_{Cr} = 52.00$ )

### Supporting electrolyte

Solution of 0.05 M DTPA (diethylenetriamino pentaacetic acid), 0.2 M  $CH_3COONa$ , 2.5 M  $NaNO_3$  at pH 6.2.

Dissolve 1.96 g of DTPA, 1.64 g of  $CH_3COONa$  and 21.3 g of  $NaNO_3$  in 80 ml of distilled water, adjust the pH to 6.2 using 30%  $NaOH$  solution.

Bring to volume in a 100ml volumetric flask with distilled water.

### Procedure

#### Analysis of Cr (VI).

Add 2 ml of supporting electrolyte to 10 ml of neutralised sample. Deaerate for 10 minutes and start the analysis.

#### Analysis of total Cr

Deaerate 10 ml of neutralised sample for 10 minutes, then add 2 ml of deaerated supporting electrolyte. Start the analysis immediately.

#### Treatment for the removal of low level of organic substances in natural water

Add 20  $\mu$ l of 30%  $H_2O_2$  to 10 ml of neutral sample irradiate for 1 hour with a 150 W UV lamp, in a closed system to avoid sample evaporation.

Then add 0.1 ml of 0.5 M  $NaOH$  and 0.5 ml of bromine water diluted 1+100. Heat at 60°C for 15 minutes.

Cool and add 0.1 ml of 0.5M  $H_2SO_4$ .

Finally, add 2 ml of supporting electrolyte.

#### Waste water, biological samples and solid samples

Opportunely digest or dissolve the sample with an oxidising treatment and dissolve the residue with distilled water. Neutralise, bring to a volume of 10 ml with distilled water and add 2 ml of supporting electrolyte.

#### Working standard solution (0.1 mg/l)

Dilute 0.1 ml of Cr (VI) concentrated standard solution in 1 l of distilled water, in a volumetric flask. Prepare this solution at the moment of the analysis.

## Warnings

Also Cr (III) give a peak in the analysis conditions, but its height decrease in few minutes, after the addition of the supporting electrolyte.

Surfactants and organic substances reduce the peak heights, than the latter must be removed using a treatment with  $H_2O_2$ . In this way only total Cr can be analysed.

The oxidation with  $H_2O_2$  has to be made at neutral pH to avoid loss of volatile forms of chromium.

Bromine water is used for the oxidation of Cr(III) to Cr(VI) because this substance give no interference with the analysis.

## Analytical report

Analysis: tap water - Chromium VI

Sample Concentration = 0.28  $\mu\text{g/l}$

Method: 3 additions

### Volumes Table

|                 |                         |
|-----------------|-------------------------|
| Solvent Volume  | 0 (ml)                  |
| Supporting Sol. | 2 (ml)                  |
| Sample Volume   | 10 (ml)                 |
| Standard Conc.  | 100 ( $\mu\text{g/l}$ ) |

### Heights Table

| # | Peak Pot. | Height              |
|---|-----------|---------------------|
| 0 | -1208.5   | 1.675 $\mu\text{A}$ |
| 1 | -1219.5   | 4.683 $\mu\text{A}$ |
| 2 | -1230.5   | 8.141 $\mu\text{A}$ |
| 3 | -1245.6   | 10.63 $\mu\text{A}$ |

### Regression Data

| # | Add. Conc.        | Height x dilution   |
|---|-------------------|---------------------|
| 0 | 0 $\mu\text{g/l}$ | 2.010 $\mu\text{A}$ |
| 1 | .5                | 5.643 $\mu\text{A}$ |
| 2 | 1                 | 9.852 $\mu\text{A}$ |
| 3 | 1.50 "            | 12.92 $\mu\text{A}$ |

$$y = ax + b$$

$$a = 7.390 \mu\text{A}^*/\mu\text{g}$$

$$b = 2.065 \mu\text{A}$$

$$r^2 = .9967$$

