

## Nickel

**Metodo alla Dimetilgliossima in tampone tartrato 0.1 M a pH 9**

**Function: Differential Pulse Voltammetry (DPV/a)**

Start Potential (mV)	-700
End Potential (mV)	-1300
Current range	2,048
Scan Speed (mV/s)	50
Number of cycles	3
Delay before sweep (s)	5
Purge and stir time (s)	300
Stirring speed (rpm)	300
Drop Size (a.u.)	60

### Nickel concentrated standard solution (1 g/l)

Dissolve 1 g of Cobalt in a minimum volume of 8 M HNO<sub>3</sub>. Bring to volume in a 1 l volumetric flask with 1% HNO<sub>3</sub>.

### Supporting Electrolyte

#### 1 – Tartrate buffer 1 M, pH 9

Dissolve 15 g of tartaric acid in 50 ml of distilled water. Add 26% NH<sub>3</sub> until pH 9.

Bring to volume with distilled water, in a 100 ml volumetric flask

#### 2 – 1% dimetilglioxime solution

Dissolve 100 mg of dimetilglioxime in 10 ml of ethanol. Prepare the solution at the moment of analysis.

#### 3- 5 M NaNO<sub>2</sub> solution

Dissolve 34.5 g of NaNO<sub>2</sub> in 100 ml of distilled water.

### Procedure

To 10 ml of sample solution add 1 ml of tartrate buffer solution, 100 µl of dimetilglioxime solution and 500 µl of NaNO<sub>2</sub> solution.

### Working standard solution (1 mg/l)

Dilute 1+999 the concentrated standard solution with distilled water. Prepare the solution at the moment of the analysis.

### Warnings

Alternative buffer solution, same pH: PIPES o HEPES.

If sample concentration is below 1 µg/l it is better to make a stripping voltammetry setting a deposition time of 30 – 120 s and a deposition potential of -700 mV.

## Analytical Report


Analysis: Nickel in deep water  
 Sample Concentration = 4.65 µg/l  
 Method: 5 additions

### Volumes Table

Solvent Volume	0 (ml)
Supporting Sol.	1.05 (ml)
Sample Volume	10 (ml)
Standard Conc.	1000 (µg/l)

### Peak Height

#	Peak Pot.	Height
0	-920.8	359.3 nA
1	-918.6	1.151 µA
2	-920.8	1.906 µA
3	-921	2.730 µA
4	-920.8	3.436 µA
5	-918.6	4.123 µA



### Regression Data

#	Add Conc.	Height x dilution	
0	0 µg/l	395.3 nA	$y = ax + b$
1	10.0 "	1.278 µA	$a = 87.43 \text{ nA}^*/\mu\text{g}$
2	20.0 "	2.136 µA	$b = 406.5 \text{ nA}$
3	30.0 "	3.085 µA	$r^2 = .9996$
4	40.0 "	3.918 µA	
5	50.0 "	4.742 µA	

