

Chloride

Function: Differential Pulse Voltammetry (DPV/a)

Start Potential (mV)	0
End Potential (mV)	400
Current range	20,048
Scan Speed (mV/s)	20
Number of cycles	3
Delay before sweep (s)	5
Purge and stir time (s)	300
Stirring speed (rpm)	300
Drop Size (a.u.)	60

Chloride concentrated standard solution (1 g/l)

Dissolve 1.6485 g of pure NaCl (dried for 2 hours at 150°C) in 1 l of distilled water, in a volumetric flask. ($MM_{\text{NaCl}} = 58.443$ $MM_{\text{Cl}} = 35.453$)

Supporting electrolyte

0.1 M KNO₃ solution

Dissolve 10.1 g of KNO₃ in 1 l of distilled water.

Procedure

Add to 10 ml of supporting electrolyte a volume of sample in way to obtain a concentration of 0.5 – 1 mg/l of Cl⁻ in the solution.

Working standard solution (100 mg/l)

Dilute 1+9 the concentrated standard solution at the moment of the analysis.

Warnings

It is not mandatory to change the internal solution of the reference electrode using the above procedure.

It is necessary to subtract the blank curve before the calculation of the peak height (point to point blank subtraction function).

Interference

Bromide and iodide have to be absent.

Analytical Report

Analysis: Tap water

Sample Concentration = 23.4 mg/l

Method: 5 additions

Blank: point to point subtraction

Volumes Table

Solvent Volume	0 (ml)
Supporting Sol.	10 (ml)
Sample Volume	0.5 (ml)
Standard Conc.	100 (mg/l)

Height Table

#	Peak Pot.	Height
0	326.1	7.205 μA
1	319.5	12.43 μA
2	314.1	18.37 μA
3	309	23.66 μA
4	304.5	28.94 μA
5	300.6	35.15 μA

Regression Data

#	Add Conc.	Height x dilution
0	0 mg/l	151.3 μA
1	20.0 "	263.7 μA
2	40.0 "	393.1 μA
3	60.0 "	511.1 μA
4	80.0 "	631.0 μA
5	100 "	773.3 μA

$$y = ax + b$$

$$a = 6.186 \mu\text{A} \cdot \text{l}/\text{mg}$$

$$b = 144.6 \mu\text{A}$$

$$r^2 = .9990$$

